

ageLOC® GALVANIC BODY SPA™



PRODUCT OVERVIEW

Your body optimized. Nu Skin has captured the future of body care with the ageLOC Galvanic Body Spa. Ergonomically designed for ease of use and handling during body treatments with ageLOC® Body Shaping Gel, the ageLOC Galvanic Body Spa with pulsating technology features a proprietary ageLOC body conductive surface. The ageLOC Galvanic Body Spa targets the sources of aging while smoothing skin to provide a more contoured, more toned appearance. It also features a self-adjusting current, lighted display, and audio indications, making it an easy to use in-home anti-aging body system.

Newly developed, patent-pending pulsating technology is optimized to maximize anti-aging effects on the arms, abdomen, buttocks, and thighs. The skin on the body is thicker and less sensitive, making it benefit from and respond to pulsating technology. It also helps stimulate, purify, and refresh your skin to reduce the visible signs of aging—for a more contoured, smoother, firmer, youthful looking you today and in the future.

TARGET AUDIENCE

PRIMARY MARKETS

- Women and men looking for a firmed, toned, younger-looking body.
- Current ageLOC Galvanic Body Spa and Galvanic Spa® customers.

SECONDARY MARKETS

- Image conscious women, ages 25–65, who may occasionally visit day spas.
- Women of any age with signs of aging on the body.

PRODUCT BENEFITS

- Optimized pulsating technology aids in stimulating skin, helps purify and refresh skin, and helps reduce the visible signs of aging.
- Laboratory research demonstrates that treatments can enhance the benefit of anti-aging ingredients for up to 24 hours.

ENHANCED FEATURES

- Automatically adjusts to your skin.
- Enhanced ergonomics—the shape of the ageLOC Galvanic Body Spa is designed to follow the contours of the body.
- ageLOC conductive surface—features a larger proprietary ageLOC conductive surface, which is able to maintain more product between the skin and the instrument surface.

USAGE

For optimal results, use your ageLOC Galvanic Body Spa with ageLOC Body Shaping Gel once a day three days a week on thighs, arms, buttocks and/or abdomen. Read through the owner's manual completely before starting your galvanic treatment.

COMPLEMENTARY PRODUCTS

- ageLOC® Body Shaping Gel—formulated to work exclusively with Nu Skin® Galvanic Spa® instruments to bring ageLOC anti-aging benefits to the body, smooth the appearance of skin, and help refresh and purify the skin while reducing the visible signs of aging.

ageLOC® GALVANIC BODY SPA™

- ageLOC® Dermatic Effects Body Contouring Lotion—scientifically formulated to bring ageLOC anti-aging benefits to the body. This daily moisturizing and contouring lotion contains ingredients that help inhibit fat production, stimulate fat breakdown, and improve the appearance while smoothing the skin.

FREQUENTLY ASKED QUESTIONS

What is unique about the ageLOC Galvanic Body Spa?

The ageLOC Galvanic Body Spa was designed for use in combination with ageLOC Body Shaping Gel on body areas—not on the face. This instrument was designed to target a larger surface area with a patent-pending ageLOC conductive surface. Secondly, this design is more ergonomic for easier use and handling during body treatments. Thirdly, the ageLOC Galvanic Body Spa and ageLOC Body Shaping Gel target end benefits/results that are different from those targeted on the face. Lastly, the skin on the body is thicker and less sensitive, making it more apt to benefit from and respond to pulsating technology.

What is the benefit of using pulsating technology on the body?

Additional research has shown that the pulsating technology enhances fluid vibration which, in turn, promotes healthy circulation. Pulsating technology aids in stimulating skin to help purify and refresh skin and help reduce the visible signs of skin aging. Pulsating technology is optimally designed for use on the body where the skin can benefit from enhanced fluid vibration.

What is the benefit of having a grooved body conductive surface?

The patent-pending grooved ageLOC conductive surface captures more key ingredients in the critical area between the conductive surface and the skin (when compared to the body conductor attachment).

How much pressure should I apply?

Nu Skin's Galvanic Spa® instruments should be very gently used with light pressure. There is no need to apply pressure to get the maximum benefits.

Should the ageLOC Galvanic Body Spa be used on the face?

No. The ageLOC Galvanic Body Spa was not designed for use on the face. This instrument was designed specifically for use on the body. It easily treats a larger surface area, easily follows the contours of the body, and, when combined with ageLOC Body Shaping Gel, brings about the targeted end results. It also provides pulsating technology, which is better suited to structural differences and targeted attributes in the area being treated. The ageLOC Galvanic Body Spa should only be used as instructed.

Are the ageLOC Galvanic Body Spa, ageLOC Body Shaping Gel, and ageLOC Dermatic Effects weight-loss products?

No, the ageLOC Galvanic Body Spa, ageLOC Body Shaping Gel, and ageLOC Dermatic Effects are not weight-loss products. ageLOC Body products were developed to smooth and promote the youthful appearance of the skin.

TO REDUCE THE RISK TO YOUR HEALTH: The ageLOC Galvanic Body Spa instrument is only for use on healthy skin. Do not use if you have open sores or wounds, skin inflammation, allergies to metal, or problematic or overly sensitive skin. Check with your physician before using if you are pregnant, have a pacemaker or similar device, are epileptic, have any metal implants, or in case of illness.

